

# Why accept bone loss?

When it comes to implant treatment, peri-implant soft tissue architecture and texture, in combination with maintained marginal bone levels, are absolute prerequisites for long lasting esthetic success. Under the right conditions, nature itself does the major part of the job. The Astra Tech implant system is designed to provide nature with the right conditions; all of the system's parts work together for reliable clinical success. In fact, the system offers unrivalled, documented results when it comes to the maintenance of marginal bone integrity and soft tissue health.

## Necessary symbiosis

The maintenance of the marginal bone is crucial both from a functional as well as an esthetic point of view. Yet, "some" bone loss is commonly accepted as an unavoidable consequence of implant treatment. Some implant suppliers even claim that it is positive and necessary in order to establish a biological width. At Astra Tech, we have never accepted this approach because there is no reason for you or your implant patients to accept bone loss. Preserving marginal bone levels and establishing the biological width at the abutment level are really all about ensuring the right stimulation of the bone and promoting

healthy soft tissue. Like the proverbial question, "Which came first, the chicken or the egg?", healthy soft tissue and maintained marginal bone are interdependent. One cannot exist without the other, because while an important task of the soft tissue is to protect the bone, the bone must be maintained to help support the soft tissue, a necessary symbiosis.


## A holistic approach

To design a successful implant system, one needs not only a great deal of knowledge about biology and mechanics, but also an understanding of what happens when the two interact. Early on in the development of the Astra Tech implant system, we realized the value of a holistic approach. That is why the Astra Tech implant system was uniquely based not only on a biological but also a biomechanical approach. With the introduction of the OsseoSpeed™ surface, we have taken this to the next level by incorporating biochemistry. Just as in nature, a successful existence cannot be determined by one single element alone. The terms biomechanics and biochemistry are no longer sufficient as there must be several interdependent features working together. This interaction is what we call the Astra Tech BioManagement Complex™.

## Astra Tech BioManagement Complex™

A successful implant system cannot be determined by one single feature alone. Just as with nature, there must be several interdependent features working together. The following combination of key features is unique to Astra Tech:

- **OsseoSpeed™** – more bone more rapidly
- **MicroThread™** – biomechanical bone stimulation
- **Conical Seal Design™** – a strong and stable fit
- **Connective Contour™** – increased soft tissue contact zone and volume



## Key references

- Hansson, S. Implant-Abutment Interface: Biomechanical Study of Flat Top versus Conical. *Clin Impl Dent Rel Res* 2000; 2: 33-41
- Norton, M. An in vitro Evaluation of the Strength of an Internal Conical Interface Compared to a Butt Joint Interface in Implant Design. *Clin Oral Impl Res* 1997; 8: 290-298
- van Steenberghe, D., De Mars G., Quirynen M., Jacobs R., Naert I. A Prospective Split-mouth Comparative Study of Two Screw-shaped Self-tapping Pure Titanium Implant Systems. *Clin Oral Impl Res* 2000; 11: 202-209
- Hansson, S. The Implant Neck: Smooth or Provided with Retention Elements. *Clin Oral Impl Res* 1999; 10: 394-405
- Hansson, S., Werke, M. The Implant Thread as a Retention Element in Cortical Bone: The Effect of Thread Size and Thread Profile: A Finite Element Study. *J Biomechanics* 2003; 36: 1247-1258
- Abrahamsson, I., Berglundh, T. Tissue Characteristics at Microthreaded Implants: An Experimental Study in Dogs. *Clin Implant Dent Rel Res* 2006; 8: 107-113
- Wennström, J. et al. Implant-supported Single-tooth Restorations: A 5-Year Prospective Study. *J Clin Periodontol* 2005; 32:567-574
- Hansson S. Surface Roughness Parameters as Predictors of Anchorage Strength in Bone: A Critical Analysis. *J Biomechanics* 2000; 33: 1297-1303
- Hansson, S., Norton, M. The Relation Between Surface Roughness and Interfacial Shear Strength for Bone-anchored Implants. A Mathematical Model. *J Biomechanics* 1999; 32: 829-836
- Gofredsen, K., Nimb, L., Hjørting-Hansen E., Jensen, J.S., Holmén, A. Histomorphometric and Removal Torque Analysis for TiO2-Blasted Titanium Implants. An Experimental Study on Dogs. *Clin Oral Impl Res* 1992; 3: 77-84
- Ivanoff, C.J., Hallgren, C., Widmark, G., Sennerby, L., Wennerberg, A. Histologic Evaluation of the Bone Integration of TiO2-Blasted and Turned Titanium Microimplants in Humans. *Clin Oral Impl Res* 2001; 12: 128-134
- Ellingsen, J. E., Johansson, C., Wennerberg, A., Holmén, A. Improved Retention and Bone-to-Implant Contact with Fluoride-Modified Titanium Implants. *Int J Oral Maxillofac Implants* 2004; 19: 659-666
- Cooper, L., Zhou, Y., Takebe, J., Guo, J., Abron, A., Holmén, A., Ellingsen, J.E. Fluoride Modification Effects on Osteoblast Behavior and Bone Formation at TiO2 Grit-blasted c.p Titanium Endosseous Implants. *Biomaterials* 2006; 27: 926 – 936
- Isa Z.M., Schneider, G., Zaharias, R., Seabold, D., Stanford, C.M. Effects of Fluoride-Modified Titanium Surfaces on Osteoblast Proliferation and Gene Expression. *Int J Oral Maxillofac Implants* 2006; 21: 203-211
- Stanford, C., Johnson, G., Fakhry, A., Gartton, D., Mellonig, J., Wagner, W. Outcomes of a Fluoride Modified Dental Implant One Year after Loading in the Posterior-Maxilla when Placed with the Osteotome Surgical Technique. *Applied Osseointegr. Res.* 2006; 5: 50 – 55
- Schliephake, H., Hüls A., Müller, M. Early Loading of Surface Modified Titanium Implants in the Posterior Mandible – Preliminary Results. *Applied Osseointegr. Res.* 2006; 5: 56 - 58
- Palmer, R.M., Palmer, P. J., Smith, B.J. A 5-Year Prospective Study of Astra Tech Single Tooth Implant. *Clin Oral Impl Res* 2000; 11: 179-182
- Puchades-Roman, L., Palmer, R.M., Palmer, P.J., Howe, L.C., Ide M., Wilson, R.F. A Clinical, Radiographic, and Microbiologic Comparison of Astra Tech and Brånemark Single Tooth Implants. *Clin Impl Dent Rel Res* 2000; 2: 78-84
- Abrahamsson, I., Berglundh, T., Wennström, J., Lindhe, J. The Peri-implant Hard and Soft Tissues at Different Implant Systems. *Clin Oral Impl Res* 1996; 7: 212-219
- Wennström, J., Ekstube, A., Grändahl, K., Karlsson, S., Lindhe, J. Implant-supported Single-tooth Restorations: A 5-Year Prospective Study. *J Clin Periodontol* 2005; 32:567-574
- De Kok, L.J., Chang, S.S., Moriarty, J.D., Cooper, L.F. A Retrospective Analysis of Peri-Implant Tissue Responses at Immediate Load/ Provisionalized Microthreaded Implants. *JOMI* 2006, 21:405-412
- Norton, M. Multiple Single-Tooth Implant Restorations in the Posterior Jaws: Maintenance of Marginal Bone Levels with Reference to the Implant-Abutment Microgap. *Int J Oral Maxillofac Implants*; in press.
- Abrahamson, I., Berglundh, T., Lindhe, J. Soft Tissue Response to Plaque Formation at Different Implant Systems. A comparative study in the Dog. *Clin Oral Impl Res* 1998;9: 73-79



$$A = \pi (r_1 \times S_1 - r_2 \times S_2)$$

$$M = F \times r$$



# Astra Tech BioManagement Complex™

—function, beauty and biology in perfect harmony





ASTRA TECH BIOMANAGEMENT COMPLEX™

### Astra Tech BioManagement Complex™

The establishment and maintenance of a soft tissue seal around the transmucosal part of an implant (i.e. the abutment) is vital for implant treatment success. The formation of the soft tissue barrier around the abutment is fundamentally a result of wound healing. Thus, during healing a barrier epithelium will form adjacent to the abutment and, apical to this epithelium, a zone of connective tissue will form and attach to the titanium surface on the abutment to protect the underlying bone tissue. The barrier epithelium and the connective tissue-implant interface will consequently establish a certain biological width of the peri-implant mucosa. It is important that this process takes place undisturbed, without micro-movements and micro-leakage in the abutment-implant connection, since this will disturb the healing process and compromise the long-term result.

Initial bone healing and long-term marginal bone stability are affected by the implant design and surface properties. Optimal biomechanical and biochemical stimuli from the implant surface are of utmost importance for the bone healing process. The long-term marginal bone stability is primarily

dependent on biomechanical stimulation from the implant, particularly around the implant neck. This means that a successful clinical result, in both a short- and long-term perspective, is related to the features of an implant. Failing to control these factors may cause problems such as black triangles between teeth and, in a worst case scenario, the implant might be lost. Infection or irritation of the soft tissue can also disturb the healing process and the long term result. These problems might be caused by factors not related to the implant as such, but to a lack of maintenance and care by the patient or by his or her general health status. Very often, the problems are caused by the same factors or circumstances that led to the initial tooth loss. Our way of safeguarding a reliable, predictable and esthetic result both in the short- and long-term, is with the Astra Tech BioManagement Complex™. This is a unique combination of the following features:


- OsseoSpeed™
- MicroThread™
- Conical Seal Design™
- Connective Contour™

The unique soft tissue seal on the abutment level, together with well-calculated biomechanical interactions with the bone around the implant neck, ensure optimal conditions for the bone. The load distribution and the lack of micro-movements and micro-leakage are the main reasons for maintained healthy peri-implant tissues and preserved marginal bone levels.



## A unique combination


To put it simply, with the Astra Tech implant system, esthetics are integrated into the implant system design. We work together with nature in supporting the natural healing process instead of interfering with it. That is why you and your patients can rely on the Astra Tech implant system, not only today, but also tomorrow and beyond.



### OsseoSpeed™ – more bone more rapidly

Building on the proven success of TiOblast™, OsseoSpeed™ is the first and only implant in the world with a chemically modified titanium


surface with a unique nanoscale topography. It stimulates early bone healing and speeds up the bone healing process. The result of the micro-roughened titanium surface treated with fluoride is increased bone formation and stronger bone-to-implant bonding. Together with MicroThread™ on the implant neck, OsseoSpeed provides true growing power in action for more reliable and effective treatment. The clinical benefits of OsseoSpeed are proven and well-documented.



### MicroThread™ – biomechanical bone stimulation

The neck of Astra Tech implants are designed with MicroThread™ — minute threads that offer optimal


load distribution and lower stress values. This design is based on a thorough understanding of bone physiology, vital to optimal implant design. Since bone tissue is designed to carry loads, dental implants must be developed to mechanically stimulate the surrounding bone in order to preserve it, taking into consideration that the critical point of the implant-bone interface is located at the marginal cortical bone where peak stresses occur.



### Conical Seal Design™ – a strong and stable fit

The Conical Seal Design™ is a conical connection below the marginal bone level that transfers the load deeper down in

the bone. Compared to conical connections above the marginal bone level and flat-to-flat designs, Conical Seal Design reduces peak stresses and thereby preserves the marginal bone. It also seals off the interior of the implant from surrounding tissues, minimizing micro-movements and micro-leakage. Conical Seal Design simplifies maintenance and ensures reliability in all clinical situations. What's more, the tight and precisely fitting implant-abutment relation of the Conical Seal Design makes abutment connection a quick and simple procedure. The abutment is self-guiding and the installation procedure is non-traumatic, eliminating the risk of bone damage.



### Connective Contour™ – increased soft tissue contact zone and volume

The Connective Contour™ is the unique contour that is created when you connect the abutment to the implant. This contour allows for an increased connective soft tissue contact zone both in height and volume, which integrates with the transmucosal part of the implant, sealing off and protecting the marginal bone.